

## Quick Coomassie Stain (QC Stain)

Cat # NB-45-00078-1L

Quick Coomassie is a new revolution in rapid 1-step Coomassie staining. The proprietary formulation, incorporating Colloidal Coomassie, is used for rapid protein staining in polyacrylamide gels.


### Storage condition

Upon receipt, store product at 4°C. Shake before use.

The QC stain is stable up to 6 months at room temperature and up to 1 year at 4°C.

### Simple 1-step Protocol

1. Pour 25 ml QC stain into a container. Use more stain if you are using a larger gel tray.
2. Remove the gel from the cassette and place the gel into the stain.
3. Leave the gel, while shaking, for a minimum of 15 minutes or until all weak protein bands are fully developed. Stain intensity is high after about 1 - 2 hours and maximum after overnight incubation.
4. Transfer the gel to DI water to remove any background staining and for gel storage. (N.B: A minimum 1 hour full stain is recommended before storing the gel in water.)



## Microwave Procedure for Gels

---

1. Using a microwave to heat up the QC stain can speed up the development of the protein bands.
2. For turbo-charging the stain, we recommend microwaving the gel, immersed in QC stain, in a suitable microwave-safe tray for a maximum 10 seconds at full power.
3. Remove the tray from the microwave and keep the gel in the QC stain for at least 30 min – 1 hour before storing the gel in DI water.

## For Mass Spectrometry Applications

---

1. Stain the gel as normal.
2. Excise the protein band of interest and put in a clean microfuge tube ideally.
3. Add 1 ml of 30% ethanol or 30% acetone.
4. Incubate for 20 min (60°C – 70°C increases the rate of de-staining).
5. Decant supernatant and repeat step 3 and 4 at least 3 times or until the gel fragment is clear.
6. Analyze the sample following the standard procedure for mass spectrometer analysis.

*For reference only*

*For Research Use Only. Not for Diagnostic or Therapeutic Use.*