

Product Components

Components	Catalog	Size-1	Size-2
		5,000 U	20,000 U
Reverse Transcriptase	RM21003	25 μ L	100 μ L
5X RT Buffer	RM21005	100 μ L	400 μ L

Product Description

ABScript Full Length Reverse Transcriptase for Single Cell is a genetically engineered M-MuLV reverse transcriptase with reduced RNase H activity and Template-switching activity. It adds adapter sequences to the 3' end of cDNA and, using Oligo (dT) VN Primer and Template-switching oligo (TSO) Primer, can initiate reverse transcription reactions from as little as 10 pg of total RNA, yielding higher full-length cDNA production. It is suitable for reverse transcription of small amounts of template and low-copy genes.

Product Source

The gene encoding the mutated M-MuLV Reverse Transcriptase (RNase H-) was expressed in Escherichia coli and purified through multiple steps.

Unit Definition

One unit (U) of activity is defined as the amount of enzyme required to incorporate 1 nmol of dTTP into acid-insoluble material in 10 minutes at 37°C, using Poly(rA)•Oligo d(T)18 as a template in a 50 μ L reaction.

Storage Temperature

-20°C

Application

- Suitable for single-cell full-length cDNA synthesis from eukaryotic mRNA or Total RNA samples containing a poly(A) tail.
- Input Sample Amount: 10 pg ~ 10 ng Total RNA (from 1~1000 lysed cells) as the starting material.

Precautions

- All buffers and enzymes should be kept on ice at all times. - Exercise caution when using various reagents to avoid cross-contamination between reagents and samples.
- Maintain a clean experimental area; wear clean gloves and masks during operations. All consumables used in the experiment, such as centrifuge tubes and pipette tips, must be RNase-free.

Instructions

1. RNA Pre-denaturation

1.1. Prepare the RNA denaturation reaction system according to the table below:

Reagent	Volume
Total RNA/cell*	10 μ L
Oligo (dT) RT Primer (10 μ M) (Self-prepared)	1 μ L
dNTP (5 mM) (Self-prepared)	1 μ L
Total volume	12 μL

1.2. Gently mix by pipetting, briefly centrifuge, and place on ice. Run the following program in a PCR machine (heated lid at 105°C):

Temperature	Time
72°C	5 min
Immediately place on ice after reaction completion*	2 min

Note: Immediately place on ice for at least 2 minutes after the reaction has concluded.

2. First Strand cDNA Synthesis

2.1. Remove the 5X RT Buffer and Reverse Transcriptase components from the -20°C freezer and place them on ice to thaw. Once the components are fully dissolved, vortex to mix thoroughly and briefly centrifuge. Prepare the first-strand cDNA synthesis reaction mixture according to the table below.

Reagent	Volumn
Products from step 1.2	12 μ L
5X RT Buffer	4 μ L
TSO Primer (20 μ M) (Self-prepared)	2 μ L
RNase Inhibitor (20 U/ μ L) (Self-prepared)	1 μ L
Reverse Transcriptase	1 μ L
Total volum	20 μL

2.2. Gently mix by pipetting up and down 15-20 times, and centrifuge briefly. Place the tube in a thermocycler, and run the following program (Lid preheated to 105°C).

Temperature	Time
42°C	90 min
70°C	10 min
12°C	Hold

2.3. After the program has completed, retrieve the first strand cDNA product. It can be directly used for dscDNA synthesis or stored at -80°C for later use.