



# PolyStain DS Kit - for Mouse and Rat antibody on Human tissue (DAB/Permanent Red)

For co-localization (Emerald/Permanent Red)

---

**NB-23-00119- 3(120 ml)**

**NB-23-00119- 2(36 ml)**

**NB-23-00119- 1(12 ml)**



**PolyStain DS Kit - for Mouse and Rat antibody on Human tissue  
(DAB/Permanent Red)**

NB-23-00119-1; NB-23-00119-2; NB-23-00119-3

**Storage: 2-8°C**

**INTENDED USE:**

The PolyStain DS Kit is designed for use with user supplied mouse and rat primary antibodies to detect two distinct antigens on human tissue or cell samples. This kit has been tested on paraffin embedded tissue. However, this kit can be used to stain frozen specimen and/or freshly prepared monolayer cell smears. Double staining is a common method used in immunohistostaining that allows detection of two distinct antigens in a single tissue. PolyStain DS Kit from NeoBiotech labs supplies the user with two polymer enzyme conjugates: antiMouse IgG (minimal cross reaction to rat) HRP polymer and anti-rat IgG (minimal cross reaction to mouse) AP polymer with two distinct substrates/chromogens, Permanent Red and DAB. Permanent Red reacts with anti-Rat AP polymer conjugate to produce the red color. DAB chromogen reacts with anti-Mouse HRP polymer conjugate to produce the brown color. A Primer step is used to increase specificity of antibody staining. PolyStain DS Kit is a nonbiotin system that avoids the extra steps involved in blocking non-specific binding due to endogenous biotin.

**KIT COMPONENTS:**

Component No.	Content	6mL Kit	36mL Kit	120mL Kit
<b>Reagent 1</b>	Rat Primer (RTU)	12mL	18mL x 2	120mL
<b>Reagent 2</b>	Rat AP Polymer (RTU)	6mL	18mL	60mL
<b>Reagent 3</b>	Mouse HRP Polymer (RTU)	6mL	18mL	60mL
<b>Reagent 4A</b>	DAB Substrate (RTU)	15mL	18mL x 2	120mL
<b>Reagent 4B</b>	DAB Chromogen (20x)	1.5mL	2mL	6mL
<b>Reagent 5A</b>	Permanent Red Substrate (RTU)	15mL	18mL x 2	120mL
<b>Reagent 5B</b>	Permanent Red Activator (5x)	3mL	7.2mL	12mL x 2
<b>Reagent 5C</b>	Permanent Red Chromogen (100x)	150µL	360µL	1.2mL
<b>Reagent 6</b>	NeoMount Universal (RTU)	15mL	18mL x 2	120mL

Gt=Goat Ms=Mouse

## RECOMMENDED PROTOCOL:

1. Fixation: To ensure the quality of the staining and obtain reproducible performance, user needs to supply appropriately fixed tissue and well prepared slides.
2. Tissue needs to be adhered to the slide tightly to avoid falling off.
3. Paraffin embedded sections must be deparaffinized with xylene and rehydrated with a graded series of ethanol before staining.
4. Cell smear samples should be made up to as much of a monolayer as possible to obtain satisfactory results.
5. Three control slides will aid the interpretation of the result: positive tissue control, reagent control (slides treated with Isotype control reagent), and negative control.
6. Proceed IHC staining: DO NOT let specimen or tissue dry from this point on.

Reagent	Staining Procedure	Incubation Time (Min.)
<b>1. Peroxidase and alkaline phosphatase Blocking</b> Supplied by user	a. Incubate slides in peroxidase and alkaline phosphatase blocking reagent (NeoPure Dual Enzyme Block NB-23-00193 was Recommended) for 10 minutes. b. Rinse the slide using distilled water at least twice.	10 min.
<b>2. HIER Pretreatment:</b> Refer to antibody data sheet.	a. Heat Induced Epitope Retrieval (HIER) may be required for primary antibody suggested by vendor. b. Wash with PBS for 2 min., 3 times.	
<b>3. Primary Antibody Mix:</b> one Mouse and one Rat primary antibody  Supplied by user	<b>Notes:</b> Investigator needs to optimize primary antibody titer and incubation time prior to double staining. a. Apply 2 drops or enough volume of Mouse and Rat primary antibody mixture to cover the tissue completely. Incubate in moist chamber for 30-60 min. Recommend 30min to shorten total protocol time. b. Wash with PBS/0.05% Tween20 for 2 minutes, 3 times.	30-60min
<b>4. Reagents 1:</b> Rat Primer (RTU)	a. Apply 2 drops or enough volume of <b>Reagent 1</b> Rat Primer to cover the tissue completely. Incubate in moist chamber for 10min. b. Wash with PBS/0.05% Tween20 for 2 minutes, 3 times	15 min.

<p><b>5. Reagents 2 &amp; 3:</b></p> <p><b>Reagents 2 :</b> Rat AP Polymer(RTU)</p> <p><b>Reagents 3:</b> Mouse HRP Polymer(RTU)</p>	<p><b>Note:</b> Make sufficient polymer mixture by adding <b>Reagent 2</b> (Rat AP Polymer) and <b>Reagent 3</b> (Mouse HRP Polymer) at 1:1 ratio, mix well.</p> <ol style="list-style-type: none"> <li>Apply 1 to 2 drops (50-100µL) of the mixture to cover each section.</li> <li>Incubate in moist chamber for 30 min.</li> <li>Wash with PBS/ 0.05% Tween20 for 2 min, 3 times.</li> <li>Rinse with distilled water.</li> </ol> <p><b>Make enough mixture for the experiment. Do not make extra volume as mixture is not stable for long term storage.</b></p>	<p>10 min</p>
<p><b>6. Reagent 4A and 4B</b></p> <p><b>Reagent 4A:</b> DAB Substrate(RTU)</p> <p><b>Reagent 4B:</b> DAB Chromogen (20x)</p>	<p><b>Note:</b> Make enough DAB mix by adding 1 drop of <b>Reagent 4B</b> (DAB Chromogen) in 1mL of <b>Reagent 4A</b> (DAB Substrate). Mix well. Use within 7 hours store at 4C.</p> <ol style="list-style-type: none"> <li>Apply 1 to 2 drops (50-100µL) of your DAB mixture to cover the tissue completely.</li> <li>Incubate for 5min.</li> <li>Rinse slides in multiple changes of distilled water 3 times, 2 each time or under running tap water for 1minute</li> </ol>	<p>5 min</p>
<p><b>7. Reagent 5A, 5B, 5C</b></p> <p><b>Reagent 5A:</b> Permanent Red Substrate (RTU)</p> <p><b>Reagent 5B:</b> Permanent Red Activator (5x)</p> <p><b>Reagent 5C:</b> Permanent Red Chromogen (100x)</p>	<ol style="list-style-type: none"> <li>Add 200µL of <b>Reagent 5B</b> (Activator) into 1mL of <b>Reagent 5A</b> (Substrate buffer) and mix well. Add 10µL of <b>Reagent 5C</b> (Chromogen) into the mixture and mix well. (<b>Note:</b> For fewer slides, Add 100µL of <b>Reagent 5B</b> (Activator) into 500µL of <b>Reagent 5A</b> (Substrate buffer) and mix well. Add 5µL of <b>Reagent 5C</b> (Chromogen) into the mixture and mix well. )</li> <li>Apply 2 drops (100µL) or enough volume of Permanent Red working solution to completely cover the tissue. Incubate for 10 min, observe appropriate color development.</li> <li>Rinse well with distilled water.</li> </ol>	<p>10 min</p>
<p><b>8. HEMATOXYLIN</b> Not provided</p>	<ol style="list-style-type: none"> <li>Counterstain with 2 drops (100µl) or enough volume of hematoxylin to completely cover tissue. Incubate for 10-15 seconds.</li> <li>Rinse thoroughly with tap water for 2-3 min.</li> <li>Put slides in PBS until show blue color (about 30 - 60 sec.)</li> <li>Rinse well in distilled water.</li> </ol>	
<p><b>9. Reagent 6:</b> NeoMount Universal (RTU)</p>	<p>Apply 2 drops (100µL) or enough volume of <b>Reagent 6</b> NeoMount Universal to cover tissue when tissue is wet. Rotate the slides to allow NeoMount Universal spread evenly.</p>	

## **PROTOCOL NOTES:**

1. The fixation, tissue slide thickness, antigen retrieval and primary antibody dilution and incubation time affect results significantly. Investigator needs to consider all factors and determine optimal conditions when interpreting the result.
2. Permanent Red is insoluble in organic solvent and can be coverslipped as well. However the dehydration steps must be shorter for optimal tissue structure and chromogen signal maintenance.
  - a. 1x 80% Ethanol 20 seconds;
  - b. 1x 95% Ethanol 20 seconds;
  - c. 3x 100% Ethanol 20 seconds each;
  - d. 1x 100% Xylene 20 seconds;
  - e. Add 1 drop of NeoMount Perm (Cat. No. Permanent mount, NB-23-00156) and coverslip.  
Press to push the air bubble out

## **PRECAUTIONS:**

DAB may be carcinogenic. Please wear gloves and take other necessary precautions

**FOR RESEARCH USE**

## Work Sheet for NB-23-00119 Kit

We designed these work sheets to help you track of each step. When staining fails these sheets help our technical support staff to pinpoint the problem. To insure that all steps are done properly, we recommend that the user fill in the actual time of their experimental step and any variation. Results will vary if time recommendations are not followed. RTU translates to ready to use.

- Used for tester to check “√” each step during the experiment
- Steps follow after de-paraffinization
- Refer to insert for details of each step

NB-23-00119 Protocol is suitable when both mouse and rat primary antibodies need or do not need pre-treatment step.

Protocol Step	NB-23-00119 Protocol	Experiment 1 Date:	Experiment 2 Date:	Experiment 3 Date:	Experiment 4 Date:
<b>Step 1</b>	Peroxidase or Alkaline Phosphatase Block User supplied				
<b>Step 2</b> Optional	HIER if needed User supplied (up to 60 min)				
<b>Step 3</b>	Mouse 1°Ab & Rat 1°Ab mix (30-60 min.)				
<b>Step 4</b>	<b>Reagent 1</b> Rat Primer RTU (10min)				
<b>Step 5</b>	<b>Reagent 2 &amp; Reagent 3</b> Rat AP Polymer & Mouse HRP Polymer require mixing (30min) Rinse with distilled water.				
<b>Step 6</b>	<b>Reagent 4A &amp; 4B</b> DAB Requires mixing! (5 min.)				

<b>Step 7</b>	<b>Reagent 5A, 5B&amp; 5C</b>  Permanent Red Requires mixing! (10 min)				
<b>Step 8</b>	Counter stain User supplied				
<b>Step 9</b>	<b>Reagent 6</b> NeoMount Universal (RTU)				
<b>Result</b>	Stain pattern on controls are correct: Fill in Yes or NO				

Result: