
Product Manual


α -Ketoglutarate Assay Kit (Fluorometric)

Catalog Number

MET-5132

100 assays

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures



CELL BIOLABS, INC.
Creating Solutions for Life Science Research

Introduction

α -Ketoglutarate (α -KG, or 2-oxoglutarate), is an important molecule with a variety of roles within multiple metabolic pathways. α -KG is produced by the oxidative decarboxylation of isocitrate via isocitrate dehydrogenase and is also the keto acid product of the deamination of glutamate via glutamate dehydrogenase. α -Ketoglutarate exerts its influence in the tricarboxylic acid (TCA) cycle as a rate determining intermediary of cellular energy metabolism. Anaplerotic reactions can refuel the cycle by generating α -Ketoglutarate from glutamate via glutamate dehydrogenase.

α -Ketoglutarate acts as a nitrogen transporter that facilitates amino acid synthesis, while also inhibiting protein degradation in muscles. L-Alanine or L-Aspartate can react with α -ketoglutarate via one of the transaminase enzymes, alanine aminotransferase (ALT) or aspartate aminotransferase (AST), to produce glutamate, which is the basis of nitrogen utilization and elimination. Amino groups from amino acids attach to α -KG where they are transported to the liver and can enter the urea cycle. These transamination and deamination pathways are also important for ammonia detoxification, as well as formation of the neurotransmitters glutamate and *gamma*-aminobutyric acid (GABA).

Other characteristics of α -Ketoglutarate include antioxidant activity by quenching of hydrogen peroxide build up and facilitating its elimination. It also generates ATP via oxidative decarboxylation for energy in intestinal cell processes and redox homeostasis. The molecule can influence hydroxylation of proline to hydroxyproline, thus increasing the formation of collagen and bone matrix formation. Recently, α -KG has been noted to extend cell lifespan and delay age-related diseases (Figure 1).




Figure 1: α -Ketoglutarate's influence on many cellular functions.

Cell Biolabs' α -Ketoglutarate Assay Kit is a simple fluorometric assay that measures the activity of α -Ketoglutarate present within plasma, serum, tissue homogenates, or cell suspensions in a 96-well microtiter plate format. Each kit provides sufficient reagents to perform up to 100 assays*, including blanks, standards and samples. Sample α -Ketoglutarate concentrations are determined by comparison with a known α -Ketoglutarate standard. The kit has a detection sensitivity limit of $\sim 0.2 \mu\text{M}$ α -Ketoglutarate.

***Note: Pyruvate present in samples can generate background signal. If pyruvate is present in your testing sample, perform 2 assays per sample reaction. Treat one sample reaction with α -KG Reaction Converter and the other without. α -Ketoglutarate is calculated from the difference in fluorescence readings from the 2 wells.**

Assay Principle

Cell Biolabs' α -Ketoglutarate Assay Kit measures α -Ketoglutarate through coupled enzymatic reactions. α -Ketoglutarate is transaminated with the production of pyruvate. Pyruvate is then detected with the fluorometric probe. Samples and standards are incubated for 60-180 minutes and then read with a standard 96-well fluorometric plate reader (Ex. 530-570 nm/Em. 590-600 nm). Samples are compared to a known concentration of α -Ketoglutarate standard within the 96-well microtiter plate format.

Related Products

1. MET-5013: Lactate Assay (Fluorometric)
2. MET-5023: Glycogen Assay (Fluorometric)
3. MET-5029: Pyruvate Assay Kit (Fluorometric)
4. MET-5055: L-Amino Acid Assay Kit (Fluorometric)
5. MET-5124: Alanine Aminotransferase (ALT) Activity assay Kit (Fluorometric)

Kit Components (shipped on blue ice)

1. α -KG Reaction Converter (Part No. 51323C): One 500 μ L tube
2. α -KG Substrate Mix (Part No. 51322B): One 750 μ L tube
3. α -Ketoglutarate Standard (Part No. 51321C): One 100 μ L tube of a 100 mM solution
4. α -KG Enzyme Mix (Part No. 51324C): One 100 μ L tube
5. Fluorescence Probe (100X) (Part No. 260003): One 100 μ L amber tube
6. HRP (Part No. 234402): One 100 μ L tube of a 100 U/mL solution in glycerol
7. 10X Assay Buffer (Part No. 50292A): One 25 mL bottle

Materials Not Supplied

1. Distilled or deionized water
2. 10 kDa MWCO centrifugal filter (e.g., Amicon Ultra 0.5 mL) (Optional) (for high protein content samples)
3. 1X PBS
4. Standard 96-well fluorescence microtiter plate

Storage

Upon receipt, store the α -KG Substrate Mix and 10X Assay Buffer at 4°C. Store the remaining components at -20°C.

Preparation of Reagents

- 1X Assay Buffer: Warm the 10X Assay Buffer to room temperature prior to using. Dilute the Assay Buffer to 1X with deionized water by diluting the 25 mL Buffer with 225 mL deionized water for 250 mL total. Mix to homogeneity. Store the 1X Assay Buffer at 4°C up to six months.

- **Reaction Reagent:** Prepare the Reaction Reagent solution by diluting the kit components accordingly. Dilute the α -KG Reaction Converter 1:20, α -KG Substrate Mix 1:20, α -KG Enzyme Mix 1:100, HRP 1:500, and Fluorescence Probe 1:100 in 1X Assay Buffer. See Table 1A below for examples of Reaction Reagent preparation based on the number of assays employed.

Note: Prepare the Reaction Reagent solution sequentially based on the chart below, beginning with the 1X Assay Buffer and adding each component to the mixture. Add the Fluorescence Probe last to the solution. Vortex the solution thoroughly between each component addition, and protect the solution from light until use. Prepare only enough for immediate use by scaling the chart examples proportionally based on the number of assays needed. For best results, place the Reaction Reagent on ice and use within 30 minutes of preparation. Do not store Reaction Reagent solutions.

1X Assay Buffer (μ L)	α -KG Reaction Converter (μ L)	α -KG Substrate Mix (μ L)	α -KG Enzyme Mix (μ L)	HRP (μ L)	Fluorescence Probe (100X) (μ L)	Number of Assays	Total Volume (100 μ L/well)
4390	250	250	50	10	50	50	5,000
2195	125	125	25	5	25	25	2,500

Table 1A. Preparation of Reaction Reagent.

Note: The Fluorescence Probe is light sensitive and must be stored accordingly.

- **Negative Control Reaction Reagent:** Prepare a Negative Control Reaction Reagent by omitting the α -KG Reaction Converter accordingly. Dilute the α -KG Substrate Mix 1:20, α -KG Enzyme Mix 1:100, HRP 1:500, and Fluorescence Probe 1:100 in 1X Assay Buffer. See Table 1B below for examples of the Negative Control Reaction Reagent preparation based on the number of assays employed.

Note: Prepare the Negative Control Reaction Reagent solution sequentially based on the chart below, beginning with the 1X Assay Buffer and adding each component to the mixture. Omit the α -KG Reaction Converter. Add the Fluorescence Probe last to the solution. Vortex the solution thoroughly between each component addition, and protect the solution from light until use. Prepare only enough for immediate use by scaling the chart examples proportionally based on the number of assays needed. For best results, place the Negative Control Reaction Reagent on ice and use within 30 minutes of preparation. Do not store Negative Control Reaction Reagent solutions.

1X Assay Buffer (μ L)	α -KG Reaction Converter (μ L)	α -KG Substrate Mix (μ L)	α -KG Enzyme Mix (μ L)	HRP (μ L)	Fluorescence Probe (100X) (μ L)	Number of Assays	Total Volume (100 μ L/well)
4640	0	250	50	10	50	50	5,000
2320	0	125	25	5	25	25	2,500

Table 1B. Preparation of Negative Control Reaction Reagent.

Note: The Fluorescence Probe is light sensitive and must be stored accordingly.

Preparation of Samples

Samples should be assayed immediately or stored at -80°C prior to performing the assay. Optimal experimental conditions for samples must be determined by the investigator. The following recommendations are only guidelines and may be altered to optimize or complement the user's experimental design. A set of serial dilutions is recommended for samples to achieve optimal assay results and minimize possible interfering compounds. Run proper controls as necessary. Always run a standard curve with samples. High protein concentrations may interfere with the assay. In this case, filter the sample with a 10kDa MWCO centrifugal filter before assaying (to reduce the possibility of protein interference and turbidity).

- **Tissues:** Weigh 500-1000 mg of sample and mince with scissors and a dounce until tissue is thoroughly liquified. Add 2 mL of 1X Assay Buffer or PBS and further sonicate the homogenate for several cycles on ice. Centrifuge 10 minutes at 12,000 x g to remove debris. Recover the supernatant and recentrifuge in a separate tube to clarify it further. Recover supernatant in a fresh eppendorf tube and incubate on ice. Prepare samples for testing and store the remaining supernatant at -80°C . Prepare further dilutions in 1X Assay Buffer.
- **Cell Suspensions:** Prepare cells at 1×10^6 cells/mL and rapidly homogenize the cell pellet with 0.2 mL cold PBS or 1X Assay Buffer. Centrifuge 10 minutes at 12,000 x g to remove debris. Recover supernatant in a fresh eppendorf tube and incubate on ice. Prepare samples for testing and store the remaining supernatant at -80°C . Prepare further dilutions in 1X Assay Buffer.
- **Serum:** Collect blood without using an anticoagulant. Allow blood to clot for 30 minutes at room temperature. Centrifuge at 2000 x g and 4°C for 10 minutes. Remove the serum layer and store on ice. Take care to avoid disturbing the white buffy layer. Aliquot samples for testing and store remaining solution at -80°C . Perform serum dilutions in 1X Assay Buffer. Perform several serial dilutions to ensure values are within the range of the standard curve.
- **Plasma:** Collect blood with heparin or citrate (EDTA could cause a quenching effect) and centrifuge at 1000 x g and 4°C for 10 minutes. Remove the plasma layer and store on ice. Take care to avoid disturbing the white buffy layer. Aliquot samples for testing and store remaining solution at -80°C . Perform plasma dilutions in 1X Assay Buffer. Perform several serial dilutions to ensure values are within the range of the standard curve.

Notes:

1. *Samples with NADH concentrations above $10 \mu\text{M}$ and glutathione concentrations above $50 \mu\text{M}$ will oxidize the Fluorescence Probe and could result in erroneous readings. To minimize this interference, it is recommended that superoxide dismutase (SOD) be added to the reaction at a final concentration of 40 U/mL.*
2. *Avoid samples containing DTT or β -mercaptoethanol since the Fluorescence Probe is not stable in the presence of thiols (above $10 \mu\text{M}$).*
3. *The Fluorescence Probe is unstable at high pH (>8.5).*

Preparation of α -Ketoglutarate Standard Curve

1. Prepare fresh α -Ketoglutarate standards by diluting the α -Ketoglutarate Standard stock from 100 mM to 1 mM in 1X PBS for a 1:100 dilution. (Example: Add 10 μL of the α -Ketoglutarate Standard stock tube to 990 μL of 1X PBS). Vortex thoroughly.
2. Continue preparing a dilution series of α -Ketoglutarate in the range of 0-25 μM by diluting the α -Ketoglutarate standard preparation according to Table 2.

Tubes	1 mM α -Ketoglutarate Standard (μ L)	1X PBS (μ L)	Final α -Ketoglutarate Concentration (μ M)	α -Ketoglutarate Quantity (nmoles/well) *
1	25	975	25	1.25
2	250 of Tube #1	250	12.5	0.625
3	250 of Tube #2	250	6.25	0.313
4	250 of Tube #3	250	3.13	0.156
5	250 of Tube #4	250	1.56	0.078
6	250 of Tube #5	250	0.78	0.039
7	250 of Tube #6	250	0.39	0.02
8	0	250	0	0

Table 2. Preparation of α -Ketoglutarate Standards.

Note: *Based on 50 μ L volume/well. Do not store diluted α -Ketoglutarate standard solutions.

Assay Protocol

Each α -Ketoglutarate standard, controls and samples should be assayed in duplicate or triplicate. A freshly prepared standard curve should be used each time the assay is performed.

****Note: Pyruvate present in samples can generate background signal. If pyruvate is present in your testing sample, perform 2 assays per sample reaction. Treat one sample reaction with α -KG Reaction Converter and the other without (1X Assay buffer will be added to the Negative Control Reaction Mix in place of α -KG Reaction Converter). α -Ketoglutarate is calculated from the difference in fluorescence readings from the 2 wells.***

1. Add 50 μ L of the diluted α -Ketoglutarate standards, controls or samples to each well of a 96-well microtiter plate.
2. Add 100 μ L of the prepared Reaction Reagent to the α -Ketoglutarate standards and to one half of the paired sample wells. Mix the well contents thoroughly.
3. Add 100 μ L of the prepared Negative Control Reaction Reagent to the remaining half of the paired sample wells, and mix the well contents thoroughly.
4. Incubate the wells for 60-180 minutes at 37°C protected from light.

Note: This assay is continuous (not terminated) and therefore may be measured at multiple time points to follow the reaction kinetics.

5. Read the fluorescence of each microwell used on a fluorescence microplate reader using 530-570 nm excitation filter and 590-600 nm emission filter.

Example of Results

The following figures demonstrate typical α -Ketoglutarate Assay Kit results. One should use the data below for reference only. This data should not be used to interpret or calculate actual sample results.




Figure 2: Example α -Ketoglutarate Standard Curve.




Figure 3: Liver Homogenates. Chicken livers were homogenized and sonicated in cold 1X Assay Buffer. Upon centrifugation, the homogenate was tested according to the assay protocol.

Calculation of Results

1. Calculate the average fluorescence values for every standard, control, and sample. Subtract the average zero standard value from itself and all standard and sample values. This is the corrected fluorescence.
2. Subtract the sample fluorescence well values without α -KG Reaction Converter ($-\alpha$ -KG RC) from the sample well values containing α -KG Reaction Converter ($+\alpha$ -KG RC) to obtain the net fluorescence difference. The fluorescence difference is due to the endogenous pyruvate activity.

$$\text{Net } \Delta\text{RFU} = (\text{RFU}_{+\alpha\text{-KG RC}}) - (\text{RFU}_{-\alpha\text{-KG RC}})$$

3. Plot the corrected fluorescence for the standards against the final concentration of the α -Ketoglutarate standards from Table 2 to determine the best curve. Use the final time point

fluorescence values to plot the α -Ketoglutarate standard curve. See Figure 2 for an example standard curve.

4. Determine the α -Ketoglutarate concentration of the samples with the equation obtained from the linear regression analysis of the standard curve. Substitute the corrected fluorescence values for each sample. Only use values within the range of the standard curve. Remember to account for dilution factors.
- 5.

$$\alpha\text{-Ketoglutarate } (\mu\text{M}) = \left[\frac{\text{Sample Corrected Fluorescence}}{\text{Slope}} \right] \times \text{Sample dilution}$$

Note: α -Ketoglutarate molecular weight is 146.11 g/mol.

References

1. Hammerqvist, F., et al. (1991) *Surgery* **109**: 28-36.
2. He, L., et al. (2015) *Curr. Protein Pept. Sci.* **16**: 576-581.
3. Hou, Y.Q., et al. (2010) *Front Biosci.* **16**: 1186-1196.
4. Tatara M., et al. (2005) *Poult. Sci.* **84(10)**: 1604-1609.
5. Yao, K., et al. (2012) *Amino Acids* **42**: 2491-2500.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

Contact Information

Cell Biolabs, Inc.
5628 Copley Drive
San Diego, CA 92111
Worldwide: +1 858 271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2019-2026: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.