

# Hematoxylin, for Automation


**Description:** Hematoxylin, for Automation is a diluted format of Mayer's hematoxylin that has been specifically formulated for Autostainers that may require a longer incubation time and a slightly lighter blue nucleus.

**Uses/Limitations:** Not to be taken internally.  
For In-Vitro Diagnostic use only.  
Histological applications.  
Do not use if reagent becomes cloudy.  
Do not use past expiration date.  
Use caution when handling reagent.  
Non-Sterile.

**Control Tissue:** Any well-fixed tissue.

**Availability:**

<u>Item #</u>	<u>Volume</u>
HAQ500	500 ml
HAQ999	1000 ml
HAQ010	10 Liters



**Storage:** Room Temperature.  
Keep container tightly closed.

**Precautions:** Avoid contact with skin and eyes.  
Wash after use.  
Observe all federal, state and local environmental regulations regarding disposal


**Manual Procedure:**


1. Re-hydrate tissue as usual.
2. Apply enough Hematoxylin to completely cover specimen and incubate for 2-10 minutes with hematoxylin according to desired intensity.
3. Apply Bluing Reagent (ScyTek catalog# BRT500) for 5-30 seconds.
4. Rinse in DI water.
5. Coverslip using media of choice.

Storage: 2° C

8° C



 ScyTek Laboratories, Inc.  
205 South 600 West  
Logan, UT 84321  
U.S.A.



EC REP  
Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem, The Netherlands