

Michel's Transport Wash Buffer

Description: This reagent is intended for use as a wash buffer for Michel's Transport Fluid. Michel's Transport Fluid is intended for use as a transport media of specimens (such as renal biopsies and lymph nodes) for immunofluorescence studies. Samples may remain in media for up to 5 days at room temperature. This reagent is not a fixative and is not suitable for transporting live cells for flow cytometry. Prior to processing, specimens should be washed in 3 changes of Michel's Transport Wash Buffer for 10 minutes each.

Uses/Limitations: Not to be taken internally.
For In-Vitro Diagnostic use only.
Histological applications.
Do not use if reagents become cloudy.
Do not use past expiration date.
Use caution when handling reagents.
Non-Sterile.



Ordering Information and Current Pricing at www.scytek.com

Availability:	Item #	Volume
	MTW125	125 ml
	MTW500	500 ml
	MTW999	1000 ml
	MTW-10000	10 Liters
	MTW-20000	20 Liters



Storage: Store at 2-8°C.

Precautions: Avoid contact with skin and eyes.
Harmful if swallowed.
Follow all Federal, State, and local regulations regarding disposal.

- Procedure:**
1. Pour Michel's Transport Wash Buffer into appropriate sized container.
 2. Remove specimen from Michel's Transport Fluid (Zeus Medium) and place in container of Michel's Transport Wash Buffer.
 3. Add enough additional Wash Buffer (if required) to completely cover specimen and allow to set for 10 minutes.
 4. Pour Wash Buffer out and refill a second time for an additional 10 minutes.
 5. Pour Wash Buffer out and refill a third time for an additional 10 minutes.
 6. Pour Wash Buffer out and proceed with immunofluorescence as usual.

- References:**
1. Vaughn Jones, S.A., et al. The use of Michel's transport medium for immunofluorescence and immunoelectron microscopy in autoimmune bullous diseases. *Journal of Cutaneous Pathology*, 1995 August; 22(4): pages 365-370.

Storage: 2°C
C

ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.

CE IVD

EC REP

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem, The Netherlands