

ASSAY DEFENDER®

Total HAMA and interference blocker

Ready-to use assay diluent for reduction of interference caused by HAMA (human anti-mouse antibodies), rheumatoid factors, nonspecific binding, cross-reactivities and matrix effects.

Interference in immunoassays can cause false positive and false negative results in routine testing. Reasons for the increasing prevalence of different kinds of interference are for example:

- population aging
- a growing number of patients with:
 - chronic autoimmune diseases
 - allergic and inflammatory diseases
- a better limit of detection in modern assays

Interferences can be classified due to their biochemical reasons and impact on assay performance.

- Interference caused by antibodies from patient samples
 e.g. HAMA (human anti-mouse antibodies), HAAA (human anti-animal antibodies),
 heterophilic antibodies and rheumatoid factors from patient samples
- Interference caused by endogenous components of the sample e.g. albumins, complement, lysozyme, fibrinogen, α -1 Antitrypsin, atypically high lipid-, salt- or sugar concentrations as well as atypical viscosities
- Interference caused by assay components

 Assay components like fluorescent or enzymatic labels can cross-react with substances from the sample or change binding properties of the assay antibodies











Fig. A: Normal sandwich assay with no interference \rightarrow correct result

Fig. B: An interfering antibody forms a bridge between capture and detection antibody → false positive result

Fig. C: Binding of the interfering antibodies to the capture antibody \rightarrow false negative result

Fig. D: Components of Assay Defender® prevent binding of the interfering antibody → correct result


Assay Defender® enables effective elimination of all kinds of interference and helps to prevent false positive or false negative results. **Distributed by:**

CliniSciences Group

Use of Assay Defender®

Assay Defender® is particularly suitable for human or animal specimens containing HAMA (human anti-mouse antibodies) or other AAA (anti-animal antibodies).

Assay Defender® is a ready-to-use dilution buffer for human or animal body fluid specimens. The previously used dilution buffer can be easily replaced by Assay Defender®. Calibrators must also be diluted in Assay Defender®. Dilution in Assay Defender® can be 1:2 or higher. A useful dilution factor for most applications is 1:10 (1 part specimen + 9 parts Assay Defender®).



Assay Defender® excels at preventing false positive results in interference samples

Serum samples from patients with rheumatoid illnesses (Sample 1 and 2) were analyzed using an established research ELISA for the detection of hAxl, a biomarker for hepatocellular carcinoma. Sample 3 is an interference-free sample from a healthy donor. Samples were diluted either in pure PBS-T/BSA, PBS-T/BSA supplemented with a commercially available HAMA Blocker (HAMA Blocker 3), commercially available ready-to-use HAMA blocking diluents (HAMA Blocker 1 and HAMA Blocker 2), or CANDOR's Assay Defender®. Assay Defender® alone can completely eliminate interferences and hence give a true negative result (< 25 ng/mL) in both interference samples. The results of interference free samples are not affected by Assay Defender®.

order numbers		
Assay Defender®	50 mL	180 050
Assay Defender®	125 mL	180 125
Assay Defender®	500 mL	180 500

Bulk quantities available on request for kit manufacturers



Distributed by:



CliniSciences Group

Austria

Company: CliniSciences GmbH Address: Sternwartestrasse 76, A-1180

Wien - Austria

Telephone: +43 720 115 580 Fax: +43 720 115 577

Email: <u>oesterreich@clinisciences.com</u> Web: <u>https://www.clinisciences.com</u>



Company: CliniSciences S.R.L Address: Avenue Stalingrad 52, 1000

Brussels - Belgium Telephone: +32 2 31 50 800 Fax: +32 2 31 50 801

Belgium

Email: belgium@clinisciences.com Web: https://www.clinisciences.com



Denmark

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064

Email: danmark@clinisciences.com Web: https://www.clinisciences.com



Finland

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064

Email: suomi@clinisciences.com Web: https://www.clinisciences.com



Company: CliniSciences S.A.S Address: 74 Rue des Suisses, 92000

Nanterre- France

Telephone: +33 9 77 40 09 09 Fax: +33 9 77 40 10 11 Email: info@clinisciences.com

Email: info@clinisciences.com Web: https://www.clinisciences.com



Germany

Company: Biotrend Chemikalien GmbH Address: Wilhelm-Mauser-Str. 41-43,

50827 Köln - Germany Telephone: +49 221 9498 320 Fax: +49 221 9498 325 Email: info@biotrend.com Web: https://www.biotrend.com



Iceland

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064

Email: island@clinisciences.com Web: https://www.clinisciences.com



Ireland

France

Company: CliniSciences Limited

Address: Ground Floor, 71 lower Baggot street Dublin D02 P593 - Ireland

Dublin D02 P593 - Ireland Telephone: +353 1 6971 146 Fax: +353 1 6971 147

Email: <u>ireland@clinisciences.com</u>
Web: <u>https://www.clinisciences.com</u>



Italy

Company: CliniSciences S.r.I Address: Via Maremmana inferiore 378 Roma 00012 Guidonia Montecelio - Italy Telephone: +39 06 94 80 56 71

Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: https://www.clinisciences.com



Netherlands

Company: CliniSciences B.V. Address: Kraijenhoffstraat 137A, 1018RG Amsterdam, Netherlands Telephone: +31 85 2082 351 Fax: +31 85 2082 353

Email: nederland@clinisciences.com Web: https://www.clinisciences.com



Company: CliniSciences AS

Address: c/o MerVerdi Munkerudtunet 10

1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: https://www.clinisciences.com



Poland

Company: CliniSciences sp.Z.o.o. Address: ul. Rotmistrza Witolda Pileckiego 67 lok. 200 - 02-781 Warszawa -Poland

Telephone: +48 22 307 0535 Fax: +48 22 307 0532

Email: polska@clinisciences.com Web: https://www.clinisciences.com



Portugal

Company: Quimigen Unipessoal LDA Address: Rua Almada Negreiros, Lote 5, Loja 14, 2615-275 Alverca Do Ribatejo - Portugal Telephone: +351 30 8808 050

Fax: +351 30 8808 052 Email: info@quimigen.com Web: https://www.quimigen.pt



Company: CliniSciences Lab Solutions Address: C/ Hermanos del Moral 13 (Bajo E), 28019, Madrid - Spain Telephone: +34 91 269 40 65 Fax: +34 91 269 40 74

Email: espana@clinisciences.com
Web: https://www.clinisciences.com



Sweden

Company: CliniSciences ApS Address: Oesterbrogade 226, st. 1, Copenhagen, 2100 - Denmark Telephone: +45 89 888 349 Fax: +45 89 884 064

Email: sverige@clinisciences.com Web: https://www.clinisciences.com



Switzerland

Company: CliniSciences AG Address: Fracht Ost Flughafen Kloten CH-8058 Zürich - Switzerland Telephone: +41 (044) 805 76 81 Fax: +41 (044) 805 76 75

Email: <u>switzerland@clinisciences.com</u>
Web: <u>https://www.clinisciences.com</u>

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511

or +44 (0) 330 684 0982 Fax: +44 (0)1753 208 899 Email: uk@clinisciences.com IWeb: https://www.clinisciences.com



Company: Biotrend Chemicals LLC Address: c/o Carr Riggs Ingram, 500 Grand Boulevard, Suite 210 Miramar Beach, FL 32550- USA

Beach, FL 32501- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: https://www.biotrend-usa.com

