

Neo CHO Feed 1, Supplement for CHO cells, w/o Insulin, w/o L-Glutamine

#Cat: NB-58-0121-100ML Size: 100ml
#Cat: NB-58-0121-500ML Size: 500ml

General Information

Neo Biotech's CHO Feed 1 and Feed 2 are feed additives especially matched to the core medium Neo CHO E. Combining those increases fed batch productivity 3- to 5-fold compared to batch culture conditions. Proposed protocols facilitate easy and efficient identification of a feeding strategy. Neo CHO Feed 1 contains concentrated amino acids, vitamins, salts, trace elements and a carbon source. Neo CHO Feed 2 is a balanced formulation of selected amino acids in a concentrated form.

Product Specifications

Appearance	Clear red solution
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping conditions.	Ambient

Formulation

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (info@neo-biotech.com) or phone (+33 9 77 40 09 09).