

## Anti-Mast Cell Tryptase antibody [AA1]

#Cat: NB-22-52467-200

Size: 0.2mg

### General Information

|                          |                                                                                                                                                                  |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                               |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Mast Cell Tryptase is suitable for use in Western Blot, Immunohistochemistry, ELISA and Immunofluorescence research applications. |
| <b>Application</b>       | WB/IHC-P/IHC-F/ELISA/IF                                                                                                                                          |
| <b>Host/Source</b>       | Mouse                                                                                                                                                            |
| <b>Reactivity</b>        | Human/Cat/Dog/Monkey/Rat                                                                                                                                         |

### Product Properties

|                            |                              |
|----------------------------|------------------------------|
| <b>Clonality</b>           | Monoclonal                   |
| <b>Clone ID</b>            | AA1                          |
| <b>Concentration</b>       | 1.0 mg/mL                    |
| <b>Conjugation</b>         | Unconjugated                 |
| <b>Purification</b>        |                              |
| <b>Dilution Range</b>      |                              |
| <b>Formulation</b>         | PBS with 0.09% sodium azide  |
| <b>Isotype</b>             | IgG1                         |
| <b>Storage Instruction</b> | Store at 2-8°C upon receipt. |

### Target Information

|                    |                                                           |
|--------------------|-----------------------------------------------------------|
| <b>Gene ID</b>     |                                                           |
| <b>Gene Symbol</b> |                                                           |
| <b>Uniprot ID</b>  |                                                           |
| <b>Immunogen</b>   | Human mast cell tryptase purified from human lung tissue. |
| <b>Immunogen</b>   |                                                           |
| <b>Region</b>      |                                                           |
| <b>Specificity</b> |                                                           |
| <b>Immunogen</b>   |                                                           |
| <b>Sequence</b>    |                                                           |

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [ RUO] licence.  
This product must not be used as for diagnostic or other medical purposes.