

Human IL-6 (Interleukin 10) Sandwich ELISA Kit

#Cat: NB-22-95536-24 Size: 24T

#Cat: NB-22-95536-48 Size: 48T

#Cat: NB-22-95536-96 Size: 96T

General Information

Product type	Kits
Short Description	This IL-6 Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids.
Applications	ELISA
Reactivity	Human

Product properties

Detection	Colorimetric
Sensitivity	0.49 pg/mL
Detection limit	0.78-50 pg/mL
Assay Time	3.5h
Storage Instruction	If unopened the kit may be stored at 2-8°C up to 1 month. If the kit will not be used within 1month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO), MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.

Target Information

Sample type	Serum, plasma and other biological fluids
Gene Symbol	IL10
Uniport ID	IL10_HUMAN

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.