

AC045

Leader in Biomolecular Solutions for Life Science



Rabbit IgG isotype control

Catalog No.: AC045

Basic Information

Observed MW

16kDa

Calculated MW

Category

SMab Recombinant Monoclonal
Antibody

Applications

IHC-P,IF/ICC,IP,ChIP,FC

Cross-Reactivity

CloneNo number

ARC5105-02

Background

Recommended Dilutions

IHC-P	1:50 - 1:200
IF/ICC	1:50 - 1:200
IP	1:50 - 1:200
FC	1:50 - 1:200
ChIP	1:50 - 1:200

Immunogen Information

Gene ID	Swiss Prot
---------	------------

Immunogen

None

Synonyms

Contact

 www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Protein A/G purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.