

Goat Anti-Human IgG+IgM+IgA antibody {Alexa Fluor 488}

#Cat : NB-22-39149-100 Size : 100µl

#Cat : NB-22-39149-1 Size : 1mL

General Information

Product Type	Secondary antibodies
Short Description	Alexa Fluor 488-conjugated goat polyclonal anti-Human IgG+IgM+IgA secondary antibody. For use in most research applications.
Applications	ELISA/IF/FC
Host/Source	Goat
Reactivity	Human

Product Properties

Clonality Clone ID	Polyclonal
Concentration	1 mg/mL
Conjugation	Alexa Fluor 488
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:100-1:1000 ELISA
Formulation	Liquid in 0.01M PBS pH7.2, 1% BSA, 50% Glycerol and 0.05% Sodium Azide
Isotype	IgG, IgM, IgA
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

Target Information

Gene ID
Gene Symbol
Uniport ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence

Alexa Fluor 350	346/442	Blue
Alexa Fluor 405	401/421	Blue
Alexa Fluor 488	496/519	Green
Alexa Fluor 532	532/553	Yellow
Alexa Fluor 555	555/565	Yellow
Alexa Fluor 568	578/603	Red/Orange
Alexa Fluor 594	590/617	Red/Orange
Alexa Fluor 633	632/647	Red
Alexa Fluor 647	650/665	Red
Alexa Fluor 660	663/690	Near IR
Alexa Fluor 680	679/702	Near IR
Alexa Fluor 750	749/775	Near IR
Alexa Fluor 790	784/814	Near IR

To use the Alexa Fluors with fluorescent imagers, use a spectral line of the blue laser diode for Alexa Fluors 405, a cyan (488nm) laser for Alexa Fluors 488, a yellow (526nm) laser for Alexa Fluor 550 or 594, and a red (633nm) laser for Alexa Fluor 649. The Alexa Fluor 680 and 790 fluors are compatible with laser-and filter- based infrared imaging instruments that emit in the 700nm, and 800nm.