

Product components

Components	Component number	Size-1 1 mL	Size-2 5 mL
dNTPs (10 Mm each)	RM20120	1 mL	5×1 mL

Product Description

dNTPs (Deoxynucleotide Solution Mix) is an equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP. The solution contains 10 mM each of dATP, dCTP, dGTP and dTTP.

It is supplied in ultrapure water as a sodium salt at pH 7.5.

Storage

-20°C

Note

Can be diluted using sterile distilled water, preferably Milli-Q® water, or sterile TE (10mM Tris-HCl, 1mM EDTA , pH 7.5).