



NIPRO Cell Culture Bag

For more information : http://www.funakoshi.co.jp/exports_contents/80156

NIPRO's cell culture bag has excellent gas permeability and is suitable for the culture of suspension cells such as lymphocytes and hybridomas.

Features

- High cell culture ability : Enabling mass cell culture with high density ($1\sim3 \times 10^6$ cells/mL) in the bag.
- High safety : Substances that may induce bad effect for cell culture, such as plasticizing agent, leaching substances and eluted substance are not contained in the bag. Each bag has been tested for leaking and pyrogen (<0.5 EU/mL).
- Sterilized by gamma irradiation.
- High transparency film : Best for the microscopic observation.
- Needleless support : Reagents and sampling can be done without needle.
- *USE LOCK TYPE SYRINGE.
- Two bags can be connected by using the connecting tube (option).
- Two sizes are available:

#87-352 Culture Bag A-350NL = 350 mL

#87-302 Culture Bag A-1000NL = 1000 mL



* NIPRO's culture bag does not contain media.

Needleless support



connecting tube (option)

Product Information

[Manufacturer : KSK]

Product Name	Size	Code	Storage
Culture Bag A-350NL (350 mL)	5 sets	87-352	RT
Culture Bag A-1000NL (1000 mL)	5 sets	87-302	RT
Connecting tube for culture bag	50 pieces	87-360	RT