

IL13 recombinant protein

#cat : NB-21-25551

Size : 100ug

OVERVIEW

PRODUCT REFERENCE	NB-21-25551
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Pro36-Asn146
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	14.54 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P35225
ALIAS/SYNONYMS	NC30, Interleukin-13

WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.