

Taq DNA Polymerase

#Cat : NB-21-25376 Size :1000U

Overview

Product Reference	NB-21-25376
Origin Species	
Fragment	
Tag	
Delivery Condition	Dry Ice

Product Information

Expression System	Prokaryotic expression
Host	Escherichia coli (E.coli)
Molecular Weight	94kDa
Form	Frozen Liquid or Lyophilized
Buffer	20mM Tris-HCl (pH8.0)、0.1 mM EDTA、1mM DTT、100 mM KCl、Stabilizers、50% glycerol
Stability & Storage	Store at -20 ° C, 12 months shelf life
Uniprot	
Alias/Synonyms	

Warning: This Product Is For Research Use Only.