

Cynomolgus Monkey Plasma, Lithium Heparin, Female individual, not filtered

#Cat: NB-53-0014-2

Size: 2ml

Lot Number(s): CYN352272-CYN352281
Expiration Date: SEE BELOW TABLE
Storage Conditions: -20°C OR COLDER
Anticoagulant: LiHep

PROJECT DETAILS

Total Volume Shipped:	20 ML as 10 x 2ML
Presentation:	CYNOMOLGUS MONKEY (EU) PLASMA LIHEP FROM FEMALE ANIMAL
Remarks:	Samples originate from Mauritian Cynomolgus macaques that are naturally B.virus free. All animals going into colonies are tested free of Herpes B. virus, SRV, SIV and STLV. All colonies are measles vaccinated. Mauritius is free of malaria and rabies. NO FREEZE/THAW CYCLES
Strain:	MACACA FASCICULARIS
Filtration:	NONE

LOT	ANIMAL UID	COLLECTION DATE	EXPERATION DATE
CYN352272	AK899C	5.NOV 21	5.NOV 26
CYN352281	VL27	14.Apr 23	14.Apr 28