
Mouse GnRH (Gonadotropin Releasing Hormone) Sandwich ELISA Kit

#Cat: NB-22-65093-48 Size: 48T
#Cat: NB-22-65093-96 Size: 96T

General Information


Product Type Kits
Short Description This GnRH Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids..
Application ELISA
Reactivity Mouse

Product Properties

Detection 15.6-1000pg/mL
Sensitivity 4.33pg/mL
Detection Limit 15.6-1000pg/mL
Assay Time 4.5 hrs
Storage Instructions Store the unopened kit in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS..

Target Information

Sample Types serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Gene Symbol [Gnrh1](#)
Uniprot ID [GON1_MOUSE](#)



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence.
This product must not be used as for diagnostic or other medical purposes.