

Product Description: CCKB antibody, human, C-term.
CCK2
Cholecystokinin type B receptor
Cholecystokinin type 1 receptor
G Protein-coupled receptor
Polyclonal antibody raised against the carboxy-terminus of the human CCKB receptor.

Product-No.: S11-9032-A-G

Immunogen: residues 428-447 coupled to KLH (NP_795344)
KLH- psiaslsrslsyttistlgpg

Form: IgG fraction supplied in PBS pH 7.4 containing 0.05% thimerosal as preservative
Antibodies were purified by peptide affinity chromatography

Quantity: 50 µg

Concentration: 0.1 mg/ml

Applications: IC, W
W = Western blot (1 µg/ml)
IC = Immunocytochemistry (1 µg/ml)

Species cross-reactivity: human

Sensitivity: 50 ng

Storage: 4 °C short term
-20 °C long term

