

## Other Animal Sera

#Cat: NB-58-0125B	Size: 100ml
#Cat: NB-58-0125A	Size: 500ml
#Cat: NB-58-0129B	Size: 100ml
#Cat: NB-58-0129A	Size: 500ml
#Cat: NB-58-0092	Size: 100ml
#Cat: NB-58-0100	Size: 100ml
#Cat: NB-58-0099	Size: 100ml

## General Information

Other animal sera are used in various applications, e.g. autologous cell culture systems or immunohistological methods.

NeoPorcine serum is often used in cultures that proliferate mycoplasma or viruses. Neo Goat and Sheep Serum are used in immunohistochemistry and to eliminate non-specific binding reactions.

### Applications:

- Autologous cell culture systems
- Immunohistochemistry
- Immunoblotting
- Immunofluorescence
- Immunoprecipitation

### Test Parameter:

- PH
- Osmolality
- Endotoxin
- Total protein
- Hemoglobin
- Mycoplasma
- Sterility

Product	Volume	Cat. No.	Storage
Neo Goat Serum	100 ml	NB-58-0125B	Store at $\leq -15^{\circ}\text{C}$ .
	500ml	NB-58-0125A	
Neo Porcine Serum	100 ml	NB-58-0129B	Store at $\leq -15^{\circ}\text{C}$ .
	500ml	NB-58-0129A	
Neo Mouse Serum	100 ml	NB-58-0092	Store at $\leq -15^{\circ}\text{C}$ .
Neo Rat Serum	100 ml	NB-58-0100	Store at $\leq -15^{\circ}\text{C}$ .
Neo Rabbit Serum	100 ml	NB-58-0099	Store at $\leq -15^{\circ}\text{C}$ .

## **Precautions and Disclaimer**

This product is for research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

## **Help Needed?**

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email ([info@neo-biotech.com](mailto:info@neo-biotech.com)) or phone +33 9 77 40 09 09.