

Mosunetuzumab [1905409-39-3]

#Cat: NB-64-42070-1mg Size: 1mg
#Cat: NB-64-42070-5mg Size: 5mg
#Cat: NB-64-42070-10mg Size: 10mg

Chemical Properties

Cas No: 1905409-39-3
Formula:
Molecular weight: 146.27 kDa
Appearance: no data available
Storage: store at low temperature
store at -20°C

Biological Description

Description	Mosunetuzumab (BCT-4465A) is a humanized antibody based on immunoglobulin G1 that targets CD20 (B cells) and CD3 (T cells). Mosunetuzumab redirects T cells to engage and eliminate malignant B cells and can be used to study relapsed or refractory (R/R) Bcell non-Hodgkin lymphoma (B-NHL) and refractory follicular lymphoma.
Targets(IC₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0068 mL	0.0342 mL	0.0684 mL
5 mM	0.0014 mL	0.0068 mL	0.0137 mL
10 mM	0.0007 mL	0.0034 mL	0.0068 mL
50 mM	0.0001 mL	0.0007 mL	0.0014 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use