

Mouse Anti-Human IgG3 antibody {FITC}

#Cat: NB-22-31354-100ul Size: 100µL

General Information

Product Type	Secondary antibodies
Short Description	FITC-conjugated mouse monoclonal anti-Human IgG3 secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Mouse
Reactivity	Human

Product properties

Clonality Clone ID	Monoclonal
Concentration	2mg/mL
Conjugation	FITC
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	IF 1:250 – 1:2000 FCM 1:250 – 1:2000
Formulation	Liquid in PBS containing 50% Glycerol, 0,5% BSA and 0.02% Sodium Azide.
Isotype	IgG3
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

Target Information

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence