

Anti-CD73 antibody [AD2] {FITC}

#Cat: NB-22-62786-20

Size: 20ug

#Cat: NB-22-62786-100

Size: 100ug

#Cat: NB-22-62786-200

Size: 200ug

General Information

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD73 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

Product properties

Detection	Monoclonal
Clone ID	AD2
Concentration	5 Mu L/Test
Conjugation	FITC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

Target Information

Gene ID
Gene Symbol
Uniport ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence

