

---

**Cynomolgus Monkey Serum, Female individual, not Filtered**

|                     |            |
|---------------------|------------|
| #Cat: NB-53-0019-1  | Size: 1ml  |
| #Cat: NB-53-0019-2  | Size: 2ml  |
| #Cat: NB-53-0019-5  | Size: 5ml  |
| #Cat: NB-53-0019-10 | Size: 10ml |

|                            |                 |
|----------------------------|-----------------|
| <b>Lot Number(s):</b>      | ETS160          |
| <b>Expiration Date:</b>    | SEE BELOW TABLE |
| <b>Storage Conditions:</b> | -20°C OR COLDER |
| <b>Anticoagulant:</b>      | NONNE           |

**PROJECT DETAILS**

|                              |   |
|------------------------------|---|
| <b>Total Volume Shipped:</b> | 1ml, 2ml, 5ml, 10ml   |
| <b>Presentation:</b>         | CYNOMOLGUS MONKEY (EU) SERUM<br>FROM 1 FEMALE ANIMAL  |
| <b>Remarks:</b>              | Samples originate from Mauritian Cynomolgus macaques that are naturally B.virus free. All animals going into colonies are tested free of Herpes B. virus, SRV, SIV and STLV. All colonies are measles vaccinated. Mauritius is free of malaria and rabies.<br>NO FREEZE/THAW CYCLES |
| <b>Strain:</b>               | MACACA FASCICULARIS   |
| <b>Filtration:</b>           | NONE  |

| LOT    | COLLECTION DATE | EXPIRATION DATE | GENDER          | ANIMAL UID'S | NUMBER OF VIALS |
|--------|-----------------|-----------------|-----------------|--------------|-----------------|
| ETS160 | 13 JUN 2019     | 13 JUN 2023     | 1 FEMALE ANIMAL | AF388        | 1               |