

### Cynomolgus Monkey Serum, Male individual, not Filtered

#Cat: NB-53-0021-1	Size: 1ml
#Cat: NB-53-0021-2	Size: 2ml
#Cat: NB-53-0021-5	Size: 5ml
#Cat: NB-53-0021-10	Size: 10ml

#### Storage conditions

-20°C OR COLDER

#### Other information

Lot Number(s):	CYN369621 [ID: AH837]- CYN369625 [ID: BY374]	
Collection Date:	09AUG2022	09AUG2027
Expiration Date:	04JUN2027	01AUG2027

PROJECT DETAILS	
Total Volume Shipped:	10 ML
Presentation:	5 LOTS OF 2 ML (1 X 2 ML EACH) OF CYNOMOLGUS MONKEY (EU) SERUM (MALE)
Remarks:	NO FREEZE/THAW CYCLES Samples originate from Mauritian Cynomolgus macaques that are naturally B.virus free. All animals going into colonies are tested free of Herpes B. virus, SRV, SIV and STLV. All colonies are measles vaccinated. Mauritius is free of malaria and rabies.
Strain:	MACACA FASCICULARIS
Filtration:	NONE

BIOHAZARD INFO: THIS MATERIAL SHOULD BE HANDLED AS IF CAPABLE OF TRANSMITTING INFECTIOUS AGENTS. PLEASE USE UNIVERSAL PRECAUTIONS. THESE MONKEYS HAVE BEEN INDIVIDUALLY TESTED FOR B-VIRUS (HERPES VIRUS SIMIAE) AND FOUND TO BE NEGATIVE. THIS LOT HAS BEEN PRODUCED FROM THOSE MONKEYS ONLY. THIS PRODUCT IS BEING SOLD FOR RESEARCH AND/OR MANUFACTURING PURPOSES ONLY. THE BIOLOGICAL SAMPLES SUPPLIED BY NEO-BIOTECH, OR ANY MATERIAL ISOLATED FROM THE SAMPLES, ARE FOR IN-VITRO RESEARCH USE ONLY AND ARE NOT TO BE USED AS A SOURCE OF MATERIAL FOR CLINICAL THERAPIES. HUMAN MATERIAL MAY BE USED IN VIVO IN ANIMALS. THE USER ASSUMES ALL RESPONSIBILITY FOR ITS USAGE AND DISPOSAL, IN ACCORDANCE WITH ALL REGULATIONS.

NEO-BIOTECH WARRANTS THAT AT THE TIME OF SHIPMENT THE PRODUCT IS FREE FROM DEFECTS IN MATERIAL AND WORKMANSHIP AND CONFORMS TO SPECIFICATIONS WHICH ACCOMPANY THE PRODUCT. NEO-BIOTECH MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR PERFORMANC.