

Human C-P (C-Peptide) Sandwich ELISA Kit

#Cat: NB-22-44104-24	Size: 24T
#Cat: NB-22-44104-48	Size: 48T
#Cat: NB-22-44104-96	Size: 96T

General Information

Product Type	Kits
Short Description	This C-P Sandwich ELISA is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human serum, plasma and other biological fluids.
Application	ELISA
Reactivity	Human

Product Properties

Detection	Colormetric
Concentration	4.69pg/mL
Detection Limit	7.81~500pg/mL
Assay Time	3.5h
Storage Instructions	If unopened the kit may be stored at 2-8°C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL

Target Information

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	
Uniprot ID	