

Product Description: mu-opioid receptor (MOR, MOP), Human, Rat, Mouse
anti-human, rat, mouse MOR antibody,
rabbit antibody raised against the carboxy-terminus of the
human, rat and mouse MOR receptor.

Category: Cells and Surface Antigens

Order number: MOR-1000-RM

Host: Rb

Description: Rabbit IgG

Applications: IHC (p), IC, WB, IP
W = Western blot (1:1000)
IC = Immunocytochemistry (1:500)
IHC (p) = Immunohistochemistry of formalin-fixed, paraffin-
embedded tissues (1:100)

Application: IHC (p), IC, WB, IP


Reactivity: Hu, Ms, Rt

Quantity: 1000 µl


Immunogen: Amino acid sequence 387-398 (LENLEAETAPLP) of rat and
mouse MOR receptor coupled to KLH

Storage: 4 °C (do not freeze)

Positive Control Tissue: brain tissue



MOR immunoreactivity in rat striatum.



MOR-ir in paraffin-embedded mouse DRG.